

VisionAire™ 5

Compact Stationary Oxygen Concentrators

Lightest, Quietest, and Most Power-Efficient at 3 & 5 LPM

Many oxygen patients around the world rely on the everyday dependability of oxygen concentrators for all of their in-home oxygen needs.

Patients and providers have long established their objectives concerning in-home oxygen. Only with the VisionAire[®] series compact concentrators, including VisionAire5, are all these objectives met.

Cost-Effective, Maintenance-Free Units

You can experience the added assurance of maintenance-free oxygen concentrators, fully backed by the expert in oxygen-generating technology. This helps ensure that you maximize the potential of your equipment inventory while greatly reducing or eliminating any unnecessary down time.

The technology of these filter-free units offers AirSep's time-proven expertise in a smaller cabinet that will please patients and their physicians, and meet your business requirements. Save time and money with the VisionAire Family!

With Safety in Mind

Important safety features include a backlined flowmeter for viewing accuracy and proper setting, a locking flow, and overdraw protection. Each concentrator is equipped with a resettable circuit breaker, a dual-function compliance/hour meter, and well-defined system alarms.

FEATURES

- △ Most power-efficient 3 LPM (175 W) & 5 LPM (290 W) unit on the market
- △ Compact, lightweight cabinet design Ñ 30 lb (13.6 kg)
- △ Near silent operation with a low sound level of just 40 dBA
- △ Maintenance-free
- △ Optional oxygen monitor for the VisionAire 5
- △ Optional air outlet for nebulizer treatments for VisionAire 5



VisionAire^a 5



VisionAire 5

Oxygen Concentration* 1-5 liters per minute at 90% + 5.5/ -3%

Dimensions In (cm) 20.8 x 14.1 x 11.5
(52.8 x 35.8 x 29.2)

Weight 30 lb; shipping weight - 37 lb
(13.6 kg; shipping weight - 16.7 kg)

Electrical 220-240 V, 50 Hz, 1.5 amps, 290 watts
230 V, 60 Hz, 1.5 amps, 290 watts

Sound Level 40 dBA

Alarms Power failure, High and low pressure
High temperature, No flow, Low
concentration (with optional oxygen monitor)

Relative Humidity Up to 95% (non-condensing)

*based on an atmospheric pressure of 14.7 psia (101 kpa) at 70; f (21 ;c). operating unit outside of operational temperature range can affect performance

WARRANTY: Three year warranty.

Ordering Information

Description

VisionAire 5 with air outlet (230 V, 60 Hz unit)

Catalog Number**

AS098-709

Phone: 0345 100 0084

E-mail: support@baro2.com