

# FreeStyle® Comfort®

The perfect fit

CAIRE®  
O<sub>2</sub> is what we do.

the  
OXYGEN  
store

## Key features and benefits



### Optimal Comfort

Let your patients experience the comfort of the body-skimming, curved silhouette design. The FreeStyle Comfort portable oxygen concentrator is lightweight and compact.



### Smart O2 delivery

Unique clinical features including autoSAT®, UltraSense®, autoDOSE, and myCAIRE™ help to ensure your patients will get the oxygen they need, when they need it.



### Proper O2 Saturation

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



### Proven Quality

Robust design features, including the sieve bed technology and critical components, allow for minimum maintenance and long life. CAIRE offers one of the best warranties available, with a 3-year warranty on the device and sieve beds.

## Specifications

### FLOW RATE

Pulse Dose 1-5  
O<sub>2</sub> Output Settings:  
Setting 1-210 mL/min  
Setting 2-420 mL/min  
Setting 3-630 mL/min  
Setting 4-840 mL/min  
Setting 5-1050 mL/min

### DIMENSIONS

10.0 x 7.3 x 3.1 in  
(25.4 x 18.5 x 7.9 cm)  
8-cell battery pack  
  
11.0 x 7.3 x 3.1 in  
(27.9 x 18.5 x 7.9 cm)  
16-cell battery pack

### WEIGHT

5 lbs (2.3 kg)  
8-cell battery pack  
  
6 lbs (2.7 kg)  
16-cell battery pack

### SOUND LEVEL

39.9 dB(A) @ Setting 2

### TRIGGER SENSITIVITY

>-0.10 cm/H<sub>2</sub>O  
(adjustable)

### BATTERY DURATION

8-cell battery  
up to 8 hours  
  
16-cell battery  
up to 16 hours

### POWER

AC Power: 100-240  
VAC, 50-60Hz  
DC Power: 11-18 VDC  
(10 max amp)  
Battery: Lithium Ion

### WARRANTY

3 years  
(Includes sieve beds)  
  
1 year on battery and  
power supplies



# FreeStyle® Comfort®

Open the door to a bigger,  
more enjoyable life

**CAIRE®**  
O<sub>2</sub> is what we do.

**the OXYGEN store**

## CARRYING BAG (Holds FreeStyle Comfort only)

**Item Number**  
MI396-1



## DESKTOP BATTERY CHARGER (Battery Not Included)

**Item Numbers**  
BT036-1  
(Type B Plug)\*



BT036-2  
(Type E Plug)\*\*



## CARRY-ALL ACCESSORY BAG (Holds Unit and Accessories)

**Item Number**  
MI372-2



**Replacement Strap  
Item Number**  
MI397-1

## BATTERY PACK

**8-Cell Item Number**  
BT034-1S



**16-Cell Item Number**  
BT035-1S



## VEHICLE POWER CORD (12 VDC Outlet)

**Item Number**  
CDO41-1



## AU POWER CORD

**Item Number**  
21347183



## BACKPACK

**Item Number**  
MI459-1



## TRAVEL ADAPTER

**Item Number** 5591-SEQ



## POWER SUPPLY WITH AC INPUT CABLE

**Item Numbers**

PW036-1S  
(Type A Plug)\*

PW036-2S  
(Type C Plug)\*



## AC INPUT CABLE

**Item Numbers**  
4997-SEQ (Type A Plug)\*  
4998-SEQ (Type C Plug)\*



## POWER SUPPLY (120W 15V) (Fits All Plug Styles)

**Item Number**  
PW036-1



## GROSS PARTICLE FILTERS

**Right Filter  
Item Number**  
FI226-1

**Left Filter  
Item Number**  
FI227-1



Phone: 0345 100 0084

E-mail: [support@baro2.com](mailto:support@baro2.com)